Reflotron Glucose

REF		SYSTEM
10744948	2 x 15	Reflotron

English

Intended use

Test for the quantitative determination of glucose in blood, serum or plasma with Reflotron systems.

Summary

The body obtains glucose from glucose-supplying carbohydrates. It is the main monosaccharide in the blood with a postprandial concentration of 5 mmol/L, and is an indispensable energy-delivering substrate for cellular functions. Glucose is broken down by glycolysis. Determination of glucose is used for the diagnosis and monitoring of disorders of carbohydrate metabolism such as diabetes mellitus, idiopathic hypoglycaemia and in pancreatic islet cell tumours.

Test principle

After application to the test strip, the sample flows into the reaction zone, where, in the case of blood samples, the separation of the erythrocytes from the plasma occurs.

D-glucose is oxidized to δ-D-gluconolactone by atmospheric oxygen in the presence of glucose oxidase (GOD). The resulting hydrogen peroxide oxidizes an indicator in the presence of peroxidase (POD). The dye formed in this manner is proportional to the glucose concentration of the sample.¹

		GOD		
glucose + O ₂	→		δ-D-gluconolactone + H ₂ O ₂	
		POD		
H ₂ O ₂ + indicator	→		dye + H ₂ O	

The glucose concentration (proportional to the dye formed) is measured kinetically at a wavelength of 642 nm and 37 °C, and is displayed after about 125 seconds in mg/dL or mmol/L.

Reagents

Components per test:

GOD (Aspergillus niger) ≥ 3.2 U; POD (horseradish) ≥ 3.2 U; 3,3',5,5'-tetramethylbenzidine (indicator) 72.6 µg; buffer.

Precautions and warnings

For in vitro diagnostic use. Exercise the normal precautions required for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

Avoid any contact with the application zone of a test strip (e.g., during pipetting of sample).

Reagent handling

Test strips are ready-for-use.

Storage and stability

Store at 2-30 °C.

Do not use the test strip after the specified expiration date.

Tightly re-cap the container immediately after removing a test strip.

Specimen collection and preparation

For specimen collection and preparation only use suitable tubes or collection containers.

Capillary blood; whole blood collected in standard sample collection tubes, as well as serum prepared from this; heparinized or EDTA-blood or heparinized or EDTA-plasma.

Use fresh capillary or venous blood immediately after collection.

Use EDTA- or heparinized blood within 10 minutes of collection. If coated single-use containers or capillary pipettes are used, please observe the stability data given by the manufacturer.

Separate serum from the cellular components immediately after coagulation - but no later than 30 minutes after collection of the blood sample. Perform the glucose determination within 2 hours.

Separate EDTA and heparinized plasma from the cellular components immediately after collection of the blood sample. Perform the glucose determination within 2 hours.²

Do not use any other anticoagulants or additives.

As glycolysis takes place in all sample materials the measured glucose value may differ from the true value in the patient even if all these conditions are observed.

Sample volume: 30 µL

Materials provided

- 1 container with 15 control strips

Materials required (but not provided)

- Reflotron instrument
 - Reflotron pipette
 - Reflotron pipette tips
 - Reflotron capillary tubes

- Reflotron Precinorm U and Reflotron Precipath U
- 0.9 % NaCl
- General laboratory equipment

Assay

For optimum performance of the assay follow the directions given in this document. Refer to the appropriate operator's manual for instrument-specific instructions.

- Remove a test strip from the container. Tightly recap the container immediately after removing a test strip.
- Peel off the aluminium protective foil, taking care not to bend the test strip.
- All Reflotron tests require a sample volume of 30 µL.
- Apply the required volume of sample onto the centre of the red application zone using a pipette (e.g., Reflotron pipette) – being careful not to touch the application zone. Avoid air-bubbles.
- Open the flap or sliding cover. Within 15 seconds of applying the sample, place the test strip onto the guide, and slide it forward horizontally until it locks into place. Close the sliding cover or flap.
- The test parameter abbreviation is shown on the display, if the test strip has been correctly inserted and the magnetic code has been read. The result is displayed depending on the setting of the instrument.

Calibration

The function curve for the Reflotron Glucose assay for converting reflectance values into concentrations is defined for each lot using the Gluco-quant (hexokinase) method from Roche Diagnostics.

The parameters of the curve are automatically transferred to the instrument via the magnetic strip during testing.

The test strips are calibrated to display plasma glucose values in line with IFCC recommendation on reporting results for blood glucose.

Quality control

For quality control, use Reflotron Precinorm U and Reflotron Precipath U.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Calculation

The glucose concentration is calculated automatically from the measurements taken, and function and conversion factors read from the magnetic strip on the lower face of each test strip. The glucose concentration is displayed in mg/dL or mmol/L, depending on the way the instrument has been set to show conventional or SI units.

Conversion factor: mg/dL x 0.0555 = mmol/L.

Limitations - interference

This test cannot be used for diagnosis or exclusion of gestational diabetes mellitus.

Methyldopa, metamazole and ascorbic acid in high or toxic concentrations can lead to measurement of low glucose values. Fluoride-containing anticoagulants should not be used. In samples from patients with monoclonal gammopathies the values measured may be too high.^{3,4,2}

The following had no influence on the results in the concentration ranges tested (criterion: recovery ± 10 % of baseline): haematocrit up to 55 %, haemolysis up to 1 %, lipaemic sera, bilirubin.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

Measuring range

10.0-600 mg/dL or 0.56-33.3 mmol/L.

If the measured glucose value is above the measuring range for the Reflotron Glucose assay (indicated by > in front of the result), the test should be repeated. If the displayed > 600 mg/dL or > 33.3 mmol/L is confirmed, the serum or plasma sample can be diluted 1 + 1 with serum or plasma having a known glucose content C₀. The true glucose value C can be calculated from the measured glucose concentration C_M using the following formula: C = 2 C_M - C₀.

Expected values

Fasting glucose adults²: 60–109 mg/dL or 3.3-6.05 mmol/L.

The glucose concentration increases by about 2 mg/dL with each decade of life.

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

Specific performance data

The data for the Reflotron Glucose assay were determined in evaluation studies. The majority of the test results were within the given ranges.

Precision

Repeatability (within-run precision):

CV (coefficient of variation) 1.8 % in the normal range, 2.8 % in the pathological range; sample material: heparinized blood.

Intermediate precision (between-day precision):

CV 2.5 % in the normal range, 2.7 % in the pathological range; sample material: control sera.

Method comparison

A comparison of the Reflotron Glucose assay (y) with the hexokinase method (x) using heparinized blood gave the following correlations:

y = 0.996 x + 4.86 and 1.034 x – 4.86; n = 50; r = 0.997.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, and the Method Sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Significant additions or changes are indicated by a change bar in the margin.

--	--

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